



health

Department:
Health
PROVINCE OF KWAZULU-NATAL

Prepared by:

Initial and Surname	Designation	Signature	Date
C van der Walt	Assistant Director		18/8/21

Reviewed by Supervisor/Operational Manager:

Initial and Surname	Designation	Signature	Date
			18/8/2021

Item details

Item description (E3009595)	Collodion Adhesive for electrodiagnostic studies. Secure and stable semi-permanent adhesive for long-term and short-term electrodiagnostic studies where secure and stable electrode application is required.
Size	Tubes of 6.5-7ml
Colour	Clear to translucent yellow
Material	Specially formulated mixture of diethyl ether and methylated spirit
Packaging (unit/box)	Tubs of 20 tubes
Functionality/performance	Collodion is the adhesive of choice for long-term or any EEG/EP studies where secure and stable electrode application is required. Collodion Adhesive is a specially formulated mixture of diethyl ether and methylated spirit for application of EEG electrodes to the head.
Purpose	Electrodiagnostic studies where secure and stable electrode application is required.
Other:	Must be compatible with electroencephalographic reusable surface cup electrodes

Approved by specifications committee chairperson:

Initial and Surname	Portfolio	Signature	Date
N.P. Mthembu	scmp		18.08.2021